

Antibodies Product Datasheet

Product Name:	ILK Antibody
Catalog Number:	R12-2925
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	ILK Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ILK.
Concentration:	1 mg/ml
Gene Name:	ILK
Other Name:	ILK-1; ILK-2; 59 kDa serine/threonine-protein kinase; P59ilk
Human Swiss Prot No.:	Q13418
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ILK Antibody antibody detects endogenous levels of ILK.
Storage and Stability:	Stable at -20°C for at least 1 year.