

Antibodies Product Datasheet

Product Name:	KDM1B Antibody
Catalog Number:	C51024
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KDM1B
Concentration:	1 mg/ml
Gene Name:	KDM1B AOF1 C6orf193 LSD2
Protein Name:	Lysine-specific histone demethylase 1B (EC 1.-.-.) (Flavin-containing amine oxidase domain-containing protein 1) (Lysine-specific histone demethylase 2)
Other Name:	
Human Swiss Prot No.:	Q8NB78
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KDM1B at Human, Mouse
Storage and Stability:	Store at -20° C/1 year