

## Antibodies Product Datasheet

<b>Product Name:</b>	MKK6 (Phospho-Ser207) Antibody
<b>Catalog Number:</b>	A7154
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC, ELISA
<b>Dilutions:</b>	MKK6 (Phospho-Ser207) Antibody
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Rabbit Polyclonal
<b>Purification:</b>	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human MKK6 around the phosphorylation site of Ser207.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	MAP2K6
<b>Other Name:</b>	Dual specificity mitogen-activated protein kinase kinase 6; kinase MKK6; MAP kinase kinase 6; MAP2K6; MAPK/ERK kinase 6; MAPKK 6; MEK6; Mitogen-activated protein kinase kinase 6; MP2K6; PRKMK6; SAPKK3
<b>Human Swiss Prot No.:</b>	P52564
<b>Formulation:</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Specificity:</b>	MKK6 (Phospho-Ser207) Antibody detects endogenous levels of MKK6 only when phosphorylated at Ser207.
<b>Storage and Stability:</b>	Stable at -20°C for at least 1 year.