
Antibodies Product Datasheet

Product Name:	ApoB Monoclonal Antibody
Catalog Number:	N0036
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human ApoB expressed in E. Coli.
Concentration:	
Gene Name:	APOB
Protein Name:	Apolipoprotein B-100
Other Name:	APOB; Apolipoprotein B-100; Apo B-100
Human Swiss Prot No.:	P04114
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	ApoB Monoclonal Antibody detects endogenous levels of ApoB protein.
Storage and Stability:	-20°C/1 year