
Antibodies Product Datasheet

Product Name:	Cytochrome P450 39A1 Antibody
Catalog Number:	R12-2681
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Cytochrome P450 39A1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 39A1.
Concentration:	1 mg/ml
Gene Name:	CYP39A1
Other Name:	Cytochrome P450 39A1, 24-hydroxycholesterol 7-alpha-hydroxylase, Oxysterol 7-alpha-hydroxylase, Hcyp39a1, CP39A
Human Swiss Prot No.:	Q9NYL5
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Cytochrome P450 39A1 Antibodyantibody detects endogenous levels of Cytochrome P450 39A1.
Storage and Stability:	Stable at -20°C for at least 1 year.