

Antibodies Product Datasheet

Product Name:	M-CSF Receptor Antibody
Catalog Number:	R12-3009
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	M-CSF Receptor Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human M-CSF Receptor.
Concentration:	1 mg/ml
Gene Name:	CSF1R
Other Name:	c-fms; CD115; CSF-1-R; CSF1R; CSFMR; FMS; Fms proto-oncogene; kinase CSFR; M-CSFR; Macrophage colony stimulating factor I receptor precursor
Human Swiss Prot No.:	P07333
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	M-CSF Receptor Antibodyantibody detects endogenous levels of M-CSF Receptor.
Storage and Stability:	Stable at -20°C for at least 1 year.