
Antibodies Product Datasheet

Product Name:	Desmoyokin (Phospho-Ser5782) Antibody
Catalog Number:	A9568
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Desmoyokin (Phospho Ser5782)
Concentration:	1 mg/ml
Gene Name:	AHNAK PM227
Protein Name:	Desmoyokin (Phospho Ser5782)
Other Name:	Neuroblast differentiation-associated protein AHNAK (Desmoyokin)
Human Swiss Prot No.:	Q09666
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Desmoyokin (Phospho Ser5782)
Storage and Stability:	-20°C/1 year