
Antibodies Product Datasheet

Product Name:	Transferrin Monoclonal Antibody(3G2)
Catalog Number:	N1581
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),ELISA
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Protein
Concentration:	1 mg/ml
Gene Name:	TF
Protein Name:	Serotransferrin (Transferrin) (Beta-1 metal-binding globulin) (Siderophilin)
Other Name:	TF; Serotransferrin; Transferrin; Beta-1 metal-binding globulin; Siderophilin
Human Swiss Prot No.:	P02787
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Transferrin protein detects endogenous levels of Transferrin
Storage and Stability:	-20°C/1 year