

Antibodies Product Datasheet

Product Name:	FoxO1A (Acetyl-Lys248) Antibody
Catalog Number:	D12261
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FoxO1A (Acetyl Lys248)
Concentration:	1 mg/ml
Gene Name:	FOXO1 FKHR FOXO1A
Protein Name:	FoxO1A (Acetyl Lys248)
Other Name:	Forkhead box protein O1 (Forkhead box protein O1A) (Forkhead in rhabdomyosarcoma)
Human Swiss Prot No.:	Q12778
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat FoxO1A (Acetyl Lys248)
Storage and Stability:	-20°C/1 year