

Antibodies Product Datasheet

Product Name:	ATIC mouse mAb
Catalog Number:	N0889
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human ATIC protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	atic
Protein Name:	
Other Name:	5 aminoimidazole 4 carboxamide 1 beta D ribonucleotide transformylase/inosinicase; 5 aminoimidazole 4 carboxamide ribonucleotide formyltransferase; 5 aminoimidazole 4 carboxamide ribonucleotide formyltransferase/IMP cyclohydrolase; 5-aminoimidazole-4-carboxa
Human Swiss Prot No.:	P31939
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ATIC and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year