
Antibodies Product Datasheet

Product Name:	UBF-1 Antibody
Catalog Number:	C10403
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human UBF1. AA range:501-550
Concentration:	1 mg/ml
Gene Name:	UBTF
Protein Name:	Nucleolar transcription factor 1
Other Name:	UBTF; UBF; UBF1; Nucleolar transcription factor 1; Autoantigen NOR-90; Upstream-binding factor 1; UBF-1
Human Swiss Prot No.:	P17480
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	UBF-1 Polyclonal Antibody detects endogenous levels of UBF-1 protein.
Storage and Stability:	-20°C/1 year