

Antibodies Product Datasheet

Product Name:	Fos B (Phospho-Ser27) Antibody
Catalog Number:	A0930
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FosB around the phosphorylation site of Ser27. AA range:12-61
Concentration:	1 mg/ml
Gene Name:	FOSB
Protein Name:	Protein fosB
Other Name:	FOSB; G0S3; Protein fosB; G0/G1 switch regulatory protein 3
Human Swiss Prot No.:	P53539
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Fos B (S27) Polyclonal Antibody detects endogenous levels of Fos B protein only when phosphorylated at S27.
Storage and Stability:	-20°C/1 year