
Antibodies Product Datasheet

Product Name:	CD14 Monoclonal Antibody
Catalog Number:	N0102
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),FCM,ELISA
Dilutions:	Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CD14 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	CD14
Protein Name:	Monocyte differentiation antigen CD14
Other Name:	CD14; Monocyte differentiation antigen CD14; Myeloid cell-specific leucine-rich glycoprotein; CD antigen CD14
Human Swiss Prot No.:	P08571
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	CD14 Monoclonal Antibody detects endogenous levels of CD14 protein.
Storage and Stability:	-20°C/1 year