

## Antibodies Product Datasheet

<b>Product Name:</b>	FLT3 (Phospho-Tyr591) Antibody
<b>Catalog Number:</b>	A9521
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:1000-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized phospho peptide around human FLT3 (Tyr591)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	FLT3 CD135 FLK2 STK1
<b>Protein Name:</b>	FLT3 (Tyr591)
<b>Other Name:</b>	Receptor-type tyrosine-protein kinase FLT3 (EC 2.7.10.1) (FL cytokine receptor) (Fetal liver kinase-2) (FLK-2) (Fms-like tyrosine kinase 3) (FLT-3) (Stem cell tyrosine kinase 1) (STK-1) (CD antigen CD135)
<b>Human Swiss Prot No.:</b>	P36888
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human Mouse FLT3 (phospho-Tyr591)
<b>Storage and Stability:</b>	-20°C/1 year