

Antibodies Product Datasheet

Product Name:	PSMD11 Antibody
Catalog Number:	R12-3378
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	PSMD11 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PSMD11.
Concentration:	1 mg/ml
Gene Name:	PSMD11
Other Name:	26S proteasome non-ATPase regulatory subunit 11; 26S proteasome regulatory subunit p44.5; 26S proteasome regulatory subunit S9; PSD11
Human Swiss Prot No.:	O00231
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	PSMD11 Antibodyantibody detects endogenous levels of PSMD11.
Storage and Stability:	Stable at -20°C for at least 1 year.