Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: Caveolin-1 (Phospho-Tyr14) Antibody

Catalog Number: A7034

Reactivity: Human, Mouse, Rat

Applications: WB, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other

applications.

Modification: Phospho Specific

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human Caveolin-1 around the phosphorylation site of Tyr14. AA range:5-

54

Concentration: 1 mg/ml

Gene Name: CAV1

Protein Name: Caveolin-1

Other Name: CAV1; CAV; Caveolin-1

Human Swiss Prot No.: Q03135

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: Phospho-Caveolin-1 (Y14) Polyclonal Antibody detects endogenous

levels of Caveolin-1 protein only when phosphorylated at Y14.

Storage and Stability: -20°C/1 year