

Antibodies Product Datasheet

Product Name:	PRMT2 Antibody
Catalog Number:	C51001
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PRMT2
Concentration:	1 mg/ml
Gene Name:	PRMT2 HMT1 HRMT1L1
Protein Name:	Protein arginine N-methyltransferase 2 (EC 2.1.1.-) (Histone-arginine N-methyltransferase PRMT2) (EC 2.1.1.125)
Other Name:	
Human Swiss Prot No.:	P55345
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PRMT2 at Human
Storage and Stability:	Store at -20° C/1 year