

Antibodies Product Datasheet

Product Name:	Apaf-1-ALT Antibody
Catalog Number:	C0128
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human APAF-1-ALT. AA range:289-338
Concentration:	1 mg/ml
Gene Name:	APAF1
Protein Name:	Apoptotic protease-activating factor 1
Other Name:	APAF1-interacting protein; APIP; APIP2; CGI-29 protein; MMRP19; dJ179L10.2; likely ortholog of mouse monocyte macrophage 19
Human Swiss Prot No.:	O14727-6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Apaf-1-ALT Polyclonal Antibody detects endogenous levels of Apaf-1-ALT protein.
Storage and Stability:	-20°C/1 year