

Antibodies Product Datasheet

Product Name:	BAF170 Antibody
Catalog Number:	C11383
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SMRC2. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	SMARCC2
Protein Name:	SWI/SNF complex subunit SMARCC2
Other Name:	SMARCC2; BAF170; SWI/SNF complex subunit SMARCC2; BRG1-associated factor 170; BAF170; SWI/SNF complex 170 kDa subunit; SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily C member 2
Human Swiss Prot No.:	Q8TAQ2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BAF170 Polyclonal Antibody detects endogenous levels of BAF170 protein.
Storage and Stability:	-20°C/1 year