
Antibodies Product Datasheet

Product Name:	LTK Antibody
Catalog Number:	R12-2979
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	LTK Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LTK.
Concentration:	1 mg/ml
Gene Name:	LTK
Other Name:	Leukocyte tyrosine kinase receptor precursor; Protein tyrosine kinase-1; TYK1
Human Swiss Prot No.:	P29376
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	LTK Antibody antibody detects endogenous levels of LTK.
Storage and Stability:	Stable at -20°C for at least 1 year.