
Antibodies Product Datasheet

Product Name:	HRP*Monoclonal mouse Anti Duck IgY
Catalog Number:	SA0127
Reactivity:	DuckIgY
Conjugate:	HRP
Applications:	WB,IHC-p,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
Modification:	Not Modification Specific
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Concentration:	4-5mg/mL
Other Name:	HRP*Monoclonal mouse Anti- Duck IgY
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year