

## Antibodies Product Datasheet

<b>Product Name:</b>	E2F-3 (Acetyl-Lys168) Antibody
<b>Catalog Number:</b>	D12061
<b>Reactivity:</b>	Human:K168, Mouse:K160, Rat:K161
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB: 1:500-10000 ELISA: 1:10000
<b>Modification:</b>	Acetyl Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized acetyl-peptide from human protein at AA range: 130-200
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	E2F3 KIAA0075
<b>Protein Name:</b>	E2F transcription factor 3
<b>Other Name:</b>	E2F3 KIAA0075
<b>Human Swiss Prot No.:</b>	O00716
<b>Formulation:</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous E2F-3 (Acetyl-Lys168) protein
<b>Storage and Stability:</b>	-20°C/1 year