

Antibodies Product Datasheet

Product Name:	MEF-2 (Acetyl-Lys403) Antibody
Catalog Number:	D12314
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MEF-2 (Acetyl Lys403)
Concentration:	1 mg/ml
Gene Name:	MEF2A MEF2
Protein Name:	MEF-2 (Acetyl Lys403)
Other Name:	Myocyte-specific enhancer factor 2A (Serum response factor-like protein 1)
Human Swiss Prot No.:	Q02078
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat MEF-2 (Acetyl Lys403)
Storage and Stability:	-20°C/1 year