

Antibodies Product Datasheet

Product Name:	AGBL4 Antibody
Catalog Number:	C15302
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AGBL4. AA range:431-480
Concentration:	1 mg/ml
Gene Name:	AGBL4
Protein Name:	Cytosolic carboxypeptidase 6
Other Name:	AGBL4; CCP6; Cytosolic carboxypeptidase 6; ATP/GTP-binding protein-like 4
Human Swiss Prot No.:	Q5VU57
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AGBL4 Polyclonal Antibody detects endogenous levels of AGBL4 protein.
Storage and Stability:	-20°C/1 year