

Antibodies Product Datasheet

Product Name:	Fnk Antibody
Catalog Number:	C18598
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PLK3. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	PLK3
Protein Name:	Serine/threonine-protein kinase PLK3
Other Name:	PLK3; CNK; FNK; PRK; Serine/threonine-protein kinase PLK3; Cytokine-inducible serine/threonine-protein kinase; FGF-inducible kinase; Polo-like kinase 3; PLK-3; Proliferation-related kinase
Human Swiss Prot No.:	Q9H4B4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fnk Polyclonal Antibody detects endogenous levels of Fnk protein.
Storage and Stability:	-20°C/1 year