
Antibodies Product Datasheet

Product Name:	HDAC1 mouse Monoclonal Antibody(4E1)
Catalog Number:	N1598
Reactivity:	Human,Rat,Mouse
Applications:	IF/ICC,IHC-p
Dilutions:	IF: 1:50-200 IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of HDAC1
Concentration:	1 mg/ml
Gene Name:	HDAC1
Protein Name:	HDAC1
Other Name:	HDAC1
Human Swiss Prot No.:	Q13547
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HDAC1 protein detects endogenous levels of HDAC1
Storage and Stability:	-20°C/1 year