

Antibodies Product Datasheet

Product Name:	C3orf39 Antibody
Catalog Number:	C50537
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C3orf39
Concentration:	1 mg/ml
Gene Name:	GTDC2 AGO61 C3orf39 EOGTL
Protein Name:	Glycosyltransferase-like domain-containing protein 2 (EC 2.4.1.255) (Extracellular O-linked N-acetylglucosamine transferase-like)
Other Name:	
Human Swiss Prot No.:	Q8NAT1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C3orf39 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year