

Antibodies Product Datasheet

Product Name:	KIF20A Antibody
Catalog Number:	R12-2951
Reactivity:	Human
Applications:	WB, IHC, ELISA
Dilutions:	KIF20A Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KIF20A.
Concentration:	1 mg/ml
Gene Name:	KIF20A
Other Name:	Kinesin-like protein KIF20A; Rabkinesin-6; Rab6-interacting kinesin-like protein; GG10_2; RAB6KIFL
Human Swiss Prot No.:	O95235
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	KIF20A Antibodyantibody detects endogenous levels of KIF20A.
Storage and Stability:	Stable at -20°C for at least 1 year.