
Antibodies Product Datasheet

Product Name:	p57kip2 (PT0208) mouse mAb
Catalog Number:	N1790
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	WB 500-2000 IHC-p 1:100-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG3, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human p57kip2
Concentration:	
Gene Name:	CDKN1C KIP2
Protein Name:	Cyclin-dependent kinase inhibitor 1C (Cyclin-dependent kinase inhibitor p57) (p57Kip2)
Other Name:	
Human Swiss Prot No.:	P49918
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human p57kip2. Heat-induced epitope retrieval (HIER) Citrate buffer of pH6.0 was highly recommended as antigen repair method in paraffin section
Storage and Stability:	-20°C/1 year