

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | FITC*Polyclonal Goat Anti Pig IgG(h+l) |
| Catalog Number: | SA0243 |
| Reactivity: | PigIgG(h+l) |
| Conjugate: | FITC |
| Applications: | WB,IHC-p,IF/ICC,ELISA |
| Dilutions: | WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000 |
| Modification: | Not Modification Specific |
| Source: | Goat |
| Purification: | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Concentration: | 4-5mg/mL |
| Other Name: | FITC*Polyclonal Goat Anti- Pig IgG(h+l) |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage and Stability: | -20°C/1 year |