

Antibodies Product Datasheet

Product Name:	Otubain-1 Antibody
Catalog Number:	C19412
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Otubain-1 . at AA range: 20-100
Concentration:	1 mg/ml
Gene Name:	OTUB1
Protein Name:	Ubiquitin thioesterase OTUB1
Other Name:	OTUB1; OTB1; OTU1; HSPC263; Ubiquitin thioesterase OTUB1; Deubiquitinating enzyme OTUB1; OTU domain-containing ubiquitin aldehyde-binding protein 1; Otubain-1; hOTU1; Ubiquitin-specific-processing protease OTUB1
Human Swiss Prot No.:	Q96FW1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Otubain-1 Polyclonal Antibody detects endogenous levels of Otubain-1 protein.
Storage and Stability:	-20°C/1 year