

Antibodies Product Datasheet

Product Name:	CXXC1 Antibody
Catalog Number:	C50490
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CXXC1
Concentration:	1 mg/ml
Gene Name:	CXXC1 CFP1 CGBP PCCX1 PHF18
Protein Name:	CpG-binding protein (CXXC-type zinc finger protein 1) (PHD finger and CXXC domain-containing protein 1)
Other Name:	
Human Swiss Prot No.:	Q9P0U4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CXXC1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year