

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | DUSP13 Antibody |
| Catalog Number: | C50568 |
| Reactivity: | Human, Mouse |
| Applications: | WB |
| Dilutions: | WB 1:500-2000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal, Rabbit,IgG |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from human DUSP13 |
| Concentration: | 1 mg/ml |
| Gene Name: | DUSP13 TMDP |
| Protein Name: | Dual specificity protein phosphatase 13 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity phosphatase SKRP4) (Testis- and skeletal-muscle-specific DSP) |
| Other Name: | |
| Human Swiss Prot No.: | Q9UII6 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of DUSP13 at Human, Mouse |
| Storage and Stability: | Store at -20° C/1 year |