

Antibodies Product Datasheet

Product Name:	HDAC7A (Phospho Ser486) Antibody
Catalog Number:	C51108
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HDAC7A (Phospho Ser486)
Concentration:	1 mg/ml
Gene Name:	HDAC7 HDAC7A
Protein Name:	Histone deacetylase 7 (HD7) (EC 3.5.1.98) (Histone deacetylase 7A) (HD7a)
Other Name:	Histone deacetylase 7 (HD7) (EC 3.5.1.98) (Histone deacetylase 7A) (HD7a)
Human Swiss Prot No.:	Q8WUI4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HDAC7A (Phospho Ser486) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year