

## Antibodies Product Datasheet

<b>Product Name:</b>	Cleaved-Caspase-5 p20 (D121) Antibody
<b>Catalog Number:</b>	L0158
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Cleaved Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human Caspase 5. AA range:102-151
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CASP5
<b>Protein Name:</b>	Caspase5
<b>Other Name:</b>	CASP5; ICH3; Caspase-5; CASP-5; ICE(rel)-III; Protease ICH-3; Protease TY
<b>Human Swiss Prot No.:</b>	P51878
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Cleaved-Caspase-5 p20 (D121) Polyclonal Antibody detects endogenous levels of fragment of activated Caspase-5 p20 protein resulting from cleavage adjacent to D121.
<b>Storage and Stability:</b>	-20°C/1 year