

Antibodies Product Datasheet

Product Name:	TGF β RIII Antibody
Catalog Number:	C0339
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TGF beta Receptor III. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	TGFBR3
Protein Name:	Transforming growth factor beta receptor type 3
Other Name:	TGFBR3; Transforming growth factor beta receptor type 3; TGF-beta receptor type 3; TGFR-3; Betaglycan; Transforming growth factor beta receptor III; TGF-beta receptor type III
Human Swiss Prot No.:	Q03167
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TGF β RIII Polyclonal Antibody detects endogenous levels of TGF β RIII protein.
Storage and Stability:	-20°C/1 year