

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | Rb Antibody |
| Catalog Number: | D12183N |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from the C-terminal region of human Rb. |
| Concentration: | 1 mg/ml |
| Gene Name: | RB1 |
| Protein Name: | Retinoblastoma-associated protein |
| Other Name: | RB1; Retinoblastoma-associated protein; p105-Rb; pRb; Rb; pp110 |
| Human Swiss Prot No.: | P06400 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | Rb Polyclonal Antibody detects endogenous levels of Rb protein. |
| Storage and Stability: | -20°C/1 year |