

Antibodies Product Datasheet

Product Name:	NPY6R Antibody
Catalog Number:	C45057
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 190-270
Concentration:	1 mg/ml
Gene Name:	NPY6R NPY1RL Y2B
Protein Name:	Putative neuropeptide Y receptor type 6 (NPY6-R) (NPY Y1-like receptor) (Putative pancreatic polypeptide receptor 2) (PP2)
Other Name:	
Human Swiss Prot No.:	Q99463
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	NPY6R Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year