
Antibodies Product Datasheet

Product Name:	Olfactory receptor 56A3 Antibody
Catalog Number:	G868
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human OR56A3. AA range:266-315
Concentration:	1 mg/ml
Gene Name:	OR56A3
Protein Name:	Olfactory receptor 56A3
Other Name:	OR56A3; OR56A3P; OR56A6; Olfactory receptor 56A3; Olfactory receptor 56A6
Human Swiss Prot No.:	Q8NH54
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Olfactory receptor 56A3 Polyclonal Antibody detects endogenous levels of Olfactory receptor 56A3 protein.
Storage and Stability:	-20°C/1 year