
Antibodies Product Datasheet

Product Name:	WEE2 Antibody
Catalog Number:	R12-3681
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	WEE2 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human WEE2.
Concentration:	1 mg/ml
Gene Name:	WEE2
Other Name:	Wee1-like protein kinase 2; Wee1-like protein kinase 1B; Wee1B kinase; WEE2; WEE1B
Human Swiss Prot No.:	P0C1S8
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	WEE2 Antibodyantibody detects endogenous levels of WEE2.
Storage and Stability:	Stable at -20°C for at least 1 year.