

Antibodies Product Datasheet

Product Name:	v-SNARE Vti1p Antibody
Catalog Number:	C19503
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human v-SNARE Vti1p.
Concentration:	1 mg/ml
Gene Name:	VTI1B
Protein Name:	Vesicle transport through interaction with t-SNAREs homolog 1B
Other Name:	VTI1B; VTI1; VTI1L; VTI1L1; VTI2; Vesicle transport through interaction with t-SNAREs homolog 1B; Vesicle transport v-SNARE protein Vti1-like 1; Vti1-rp1
Human Swiss Prot No.:	Q9UEU0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	v-SNARE Vti1p Polyclonal Antibody detects endogenous levels of v-SNARE Vti1p protein.
Storage and Stability:	-20°C/1 year