

Antibodies Product Datasheet

Product Name:	FNDC5 Antibody
Catalog Number:	C51051
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FNDC5
Concentration:	1 mg/ml
Gene Name:	FNDC5 FRCP2
Protein Name:	Fibronectin type III domain-containing protein 5 (Fibronectin type III repeat-containing protein 2) [Cleaved into: Irisin]
Other Name:	
Human Swiss Prot No.:	Q8NAU1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FNDC5 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year