

Antibodies Product Datasheet

Product Name:	Anti-Human CD21 (ABT083) Mouse Monoclonal Antibody
Catalog Number:	V0229
Clone ID:	ABT083
Reactivity:	Human
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	CR2
Gene Name:	P20023
Other Name:	C3DR, CD 21, CD21, Complement C3d receptor 2, Complement C3d receptor, Complement component (3d/Epstein Barr virus) receptor 2, Complement receptor type 2, CR, Cr2, CR2_HUMAN, CVID7, EBV receptor, EBV-R, Epstein Barr virus receptor, Epstein-Barr virus receptor, EVBR, SLEB9
Isotype	IgG1, Kappa
Localization:	Tonsil/ Appendix
Positive Tissue Control:	Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CD21 protein.
Storage and Stability:	-20°C/1 year