

Antibodies Product Datasheet

Product Name:	MTMR3 Antibody
Catalog Number:	C50128
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MTMR3
Concentration:	1 mg/ml
Gene Name:	MTMR3 KIAA0371 ZFYVE10
Protein Name:	Myotubularin-related protein 3 (EC 3.1.3.48) (FYVE domain-containing dual specificity protein phosphatase 1) (FYVE-DSP1) (Zinc finger FYVE domain-containing protein 10)
Other Name:	
Human Swiss Prot No.:	Q13615
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MTMR3 at Human,Rat
Storage and Stability:	Store at -20° C/1 year