
Antibodies Product Datasheet

Product Name:	E2A Monoclonal Antibody
Catalog Number:	N0203
Reactivity:	Human
Applications:	WB,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human E2A expressed in E. Coli.
Concentration:	
Gene Name:	TCF3
Protein Name:	Transcription factor E2-alpha
Other Name:	TCF3; BHLHB21; E2A; ITF1; Transcription factor E2-alpha; Class B basic helix-loop-helix protein 21; bHLHb21; Immunoglobulin enhancer-binding factor E12/E47; Immunoglobulin transcription factor 1; Kappa-E2-binding factor; Transcription facto
Human Swiss Prot No.:	P15923
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	E2A Monoclonal Antibody detects endogenous levels of E2A protein.
Storage and Stability:	-20°C/1 year