

Antibodies Product Datasheet

Product Name:	PTEN Antibody
Catalog Number:	C30793
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 IHC-p: 100-300.WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from PTEN at AA range: 251-300
Concentration:	1 mg/ml
Gene Name:	PTEN MMAC1 TEP1
Protein Name:	PTEN
Other Name:	phosphatase and tensin homolog; phosphatase and tensin homolog pseudogene 1
Human Swiss Prot No.:	P60484
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PTEN Polyclonal Antibody detects endogenous levels of PTEN
Storage and Stability:	-20°C/1 year