
Antibodies Product Datasheet

Product Name:	CNTF Antibody
Catalog Number:	C40099
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-1000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CNTF
Concentration:	1 mg/ml
Gene Name:	CNTF
Protein Name:	Ciliary neurotrophic factor (CNTF)
Other Name:	CNTF; Ciliary neurotrophic factor; CNTF; YM3544
Human Swiss Prot No.:	P26441
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CNTF protein(A218) detects endogenous levels of CNTF
Storage and Stability:	-20°C/1 year