

Antibodies Product Datasheet

Product Name:	OLFM4 Antibody
Catalog Number:	C49766
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human OLFM4
Concentration:	1 mg/ml
Gene Name:	OLFM4 GW112 UNQ362/PRO698
Protein Name:	Olfactomedin-4 (OLM4) (Antiapoptotic protein GW112) (G-CSF-stimulated clone 1 protein) (hGC-1) (hOLF4)
Other Name:	
Human Swiss Prot No.:	Q6UX06
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of OLFM4 at Human
Storage and Stability:	Store at -20° C/1 year