

Antibodies Product Datasheet

Product Name:	SEMA3A Antibody
Catalog Number:	R12-2351
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000 ,Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SEMA3A. AA range:688-737
Concentration:	1 mg/ml
Gene Name:	SEMA3A
Protein Name:	Semaphorin-3A
Other Name:	SEMA3A; SEMAD; Semaphorin-3A; Semaphorin III; Sema III
Human Swiss Prot No.:	Q14563
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SEMA3A Polyclonal Antibody detects endogenous levels of SEMA3A protein.
Storage and Stability:	-20°C/1 year