
Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | Cytokeratin 17 Monoclonal Antibody |
| Catalog Number: | N0180 |
| Reactivity: | Human |
| Applications: | WB,IHC-p,IF(paraffin section),ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Unmodified/ Total |
| Source: | Monoclonal Mouse |
| Purification: | Affinity purification |
| Immunogen: | Purified recombinant fragment of Cytokeratin 17 expressed in E. Coli. |
| Concentration: | |
| Gene Name: | KRT17 |
| Protein Name: | Keratin type I cytoskeletal 17 |
| Other Name: | KRT17; Keratin; type I cytoskeletal 17; 39.1; Cytokeratin-17; CK-17; Keratin-17; K17 |
| Human Swiss Prot No.: | Q04695 |
| Formulation: | Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol. |
| Specificity: | Cytokeratin 17 Monoclonal Antibody detects endogenous levels of Cytokeratin 17 protein. |
| Storage and Stability: | -20°C/1 year |